

Genome Assembly Pipeline: Tools & Steps

Overview

This document outlines the tools and steps used in the **genome assembly pipeline** after data download. The steps include quality control, read trimming, assembly, evaluation, and annotation.

Table of Contents

- Overview
- Tools Required
- Pipeline Steps
 - 1. Quality Control with FastQC
 - 2. Read Trimming with Trimmomatic
 - 3. Error Correction with SPAdes
 - 4. Genome Assembly using SPAdes
 - 5. Assembly Evaluation with QUAST
 - 6. Genome Annotation with Prokka

Tools Required

The following tools must be installed before running the pipeline:

- **FastQC**: Quality control of raw and trimmed reads.
- **Trimmomatic**: Read trimming and quality filtering.
- **SPAdes**: Genome assembly and error correction.
- **QUAST**: Assembly quality assessment.
- **Prokka**: Genome annotation.

Pipeline Steps

1. Quality Control with FastQC

Command:

```
fastqc SRR9620862_1.fastq SRR9620862_2.fastq
```

- This step checks read quality, GC content, and adapter contamination.
 - Generates an HTML report for visual inspection.
-

2. Read Trimming with Trimmomatic

Command:

```
trimmomatic PE SRR9620862_1.fastq SRR9620862_2.fastq \  
SRR9620862_1_paired.fastq SRR9620862_1_unpaired.fastq \  

```

```
SRR9620862_2_paired.fastq SRR9620862_2_unpaired.fastq \
LEADING:10 TRAILING:10 SLIDINGWINDOW:5:20 MINLEN:250
```

- Trims low-quality bases and removes short reads.
 - **Paired-end reads** are handled separately to ensure both ends are retained.
-

3. Error Correction with SPAdes

Command:

```
spades.py -1 SRR9620862_1.fastq -2 SRR9620862_2.fastq \
-o spades_corrected --only-error-correction
```

- Performs read error correction before assembly.
-

4. Genome Assembly using SPAdes

Default Run

```
spades.py -1 spades_corrected/corrected/SRR9620862_100.0_0.cor.fastq.gz \
-2 spades_corrected/corrected/SRR9620862_200.0_0.cor.fastq.gz \
-o spades_default_assembly --only-assembler
```

Careful Run with K-mer Optimization

```
spades.py -k 21,33,55,77,99,127 --careful --only-assembler \
-1 spades_corrected/corrected/SRR9620862_100.0_0.cor.fastq.gz \
-2 spades_corrected/corrected/SRR9620862_200.0_0.cor.fastq.gz \
-o spades_careful_assembly
```

- The **default run** performs standard assembly.
 - The **careful run** uses optimized k-mers and error correction.
-

5. Assembly Evaluation with QUAST

Command:

```
quast -o quast_SRR9620862_out \
-R Reference_Genome/reference_genome.fna.gz \
-g Reference_Genome/anno_reference_genome.gff.gz \
--labels Spades_Default,Spades_Careful \
spades_default_assembly/contigs.fasta \
spades_careful_assembly/contigs.fasta
```

- Compares assemblies to the reference genome.

- Reports N50, GC content, number of contigs, and misassemblies.
-

6. Genome Annotation with Prokka

Command:

```
prokka --force --outdir prokka_annotation --prefix annotation \
      spades_default_assembly/contigs.fasta
```

- Identifies genes, proteins, and functional elements in the assembled genome.
 - Outputs annotation files including .gff, .gbk, and .faa.
-